

SAFETY CHEM


CHEMICAL LABELING PROGRAM

- Easily label all chemicals to meet CAP requirements
- No more hand writing each label
- Prints one label at time or full sheets as necessary
- Software includes a database of over 200 Histo Chemicals
- Easily add, edit, and save new chemicals, cautions and treatments to the database
- Perfect for Clinical or Industrial laboratories



Package includes:

Safety Chem™ Software
and set of blank
Diamond Hazard Labels



Actual size

Call us to request a complete catalogue of H/I products for HistoPathology, 1-800-442-2537



HACKER Instruments & Industries Inc.
1132 Kincaid Bridge Road
Winnsboro, SC 29180 USA
Tel: 803-712-6100 : Fax: 803-712-6116

www.hackerinstruments.com email: hackerlab@aol.com