

# HistoPro<sup>®</sup> 3030

## H&E Slide Stainer



**MADE IN THE  
USA**

*HistoPro<sup>®</sup> 3030 offers continuous load, multi-protocol, automated slide staining in a compact, price-conscious package.*

- *Compact Footprint: Only 34" (86 cm) wide*
- *30 Stations*
- *Continuous Load with Multiple Protocols*
- *5 Rinse Stations*
- *2 Water Inlets*
- *State-of-the-Art, User-Friendly Operator Interface with Touch Screen*
- *Built-In Fume Extraction Can Connect to Lab Fume Exhaust System*
- *Designed and Manufactured in the USA*

**Hacker**  
Instruments & Industries Inc.  
*Since 1942*

Visit us online @ [www.hackerinstruments.com](http://www.hackerinstruments.com)

1132 Kincaid Bridge Road • PO BOX 1176 • Winnsboro, SC 29180  
(803) 712-6100

## HistoPro 3030 Slide Stainer Specifications

### Capacity

Number of Stations:	30
Number of Rinse Stations:	5 (user selectable)
Rinse Water Flow Control:	Individually adjustable for 4 stations
Number of Water Inlets:	2
Number of Load Station:	1
Number of Unload Station:	2
Unload Station with Xylene for H & E Stains:	1
Unload Station with Water for Special Stains:	1
Reagent Container Volume:	350 ml
Slide Carrier Capacity:	30 slides
Number of Drying Stations:	2
Drying Station Temperature Range:	40° C - 70° C

### User Interface

Display:	Active matrix TFT 6.5" display
Touch Screen:	Resistive touch pad
Number of Programs:	50
Number of Steps in a Program:	100
Program Download/Upload:	Via USB port

### Physical Characteristics

Width:	34" (86.4 cm)
Depth:	25" (63.5 cm)
Height:	19" (48.3 cm)
Weight:	70 lbs (32 kg)

### Electrical Characteristics

Voltage:	115 Volts AC or 230 volts AC, 50/60 Hz
Power	1500 Watts maximum

Cat. No. 97 00 010

Hacker HistoPro® 3030 Slide Stainer

\$26,500.00

**Hacker**  
Instruments & Industries Inc.  
*Since 1942*

Visit us online @ [www.hackerinstruments.com](http://www.hackerinstruments.com)

1132 Kincaid Bridge Road • PO BOX 1176 • Winnsboro, SC 29180  
(803) 712-6100