



HistoPro 414

Linear Slide Stainer

Features

- Offers consistent staining of slides in an easy-to-use, compact, cost conscious package.
- Ideal for Mohs techniques, special stains, and general surgical pathology frozen sections.
- Fits on top of most cryostats
- Continuous loading feature to keep up with lab demand.
- Programmable station timing for 2 to 300 seconds.
- 14 Reagent stations including 3 running water stations.
- Holding station for 16 processed slides.
- Only 50 ml per reagent container.
- Up to 900 slides per hour at 10 second station time.
- Complete With: 4 Slide carriers for 4 slides each; 14 Reagent station containers (50 ml); 3 Running water rinse containers (50 ml); External water supply hose; Internal rinse station hose; Operators manual; A/C Adapter, 12V DC at 2 Amps;

<u>Cat. No.</u>		<u>Price USD</u>
414 00 01	HistoPro 414 Linear Slide Stainer	\$4,850.00
	Additional Accessories	
414 00 05	Slide Carrier	\$90.00
414 00 06	Reagent Container	\$69.00
414 00 07	Rinse Container	\$90.00